

NEW FEATURES UPDATE



V I S A R I S

TECHNICAL BULLETIN

Q2 2018

AVANSE DR VERSION 4.7.1

I-RAY DETECTOR SUPPORT

LONGER BUCKY TRAVEL IN BUCKY TABLES

**BUTTONLESS CONTROL OF WALL
STAND ELEVATION**

NEW STITCHING PLATFORM

**REAL EXPOSURE TECHNIQUE ENTRY
FOR RETRO-FIT SYSTEMS**

AVANSE DR VERSION 4.7.1

- Easier operator selection with drop menu on registration screen
 - keep track of multiple operators easily



- Auto-advance to next view in exam
 - faster operation in high throughput situations

AVANSE DR VERSION 4.7.1

- Easier APR/procedure selection on patient registration:
view only procedures or customisable ordering - loose less time on admin
- Search Modality worklists for multiple modality tags:
DX, CR, DR ... for complicated departments with multiple devices




- Clearer PACS sending status display:
quickly and efficiently

i RAY DETECTOR SUPPORT


- New detector option with good image quality at competitive price
- Available in 43x43cm fixed and two WiFi sizes (35x43cm and 43x43cm)



LONGER BUCKY TRAVEL IN BUCKY TABLES



- More than 15% longer bucky travel on 4/6 way tables
- Longer patient coverage on the table: less patient manipulation
- Longer stitching coverage on supine patients for entire spine imaging



BUTTONLESS CONTROL OF WALL STAND ELEVATION



- Intuitive buttonless brake release on wall stand elevation
- No need to hunt for release button: easier and faster positioning
- Ergonomic shape that is easy to grip in all arm elevation positions

NEW STITCHING PLATFORM

- Smoother shape, fewer corners and easy access brakes on wheels
- Lower step for easier and safer patient access
- Fold up step for spine imaging: no need for patient to climb onto the platform
- Positioning guides for manually positioned stitching



REAL EXPOSURE TECHNIQUE ENTRY FOR RETRO-FIT SYSTEMS

- Keep track of dose, even on retro-fit systems
- Enter exposure technique that is reported in DICOM

